

Instructions for Using the Percussion Pro

Go to [HumanLocomotion.com](https://www.HumanLocomotion.com) for a link to the video with full instructions for use.

Thanks for purchasing the Percussion Pro! We wanted to offer a percussion gun that was small and portable, but also offered the 50-60 cycles per second required to suppress spindle information and reduce muscle tone. Dr. Michaud also wanted a massage head attachment that was large and slightly rounded, so it would distribute pressure over a broader area and glide smoothly over the skin. When he couldn't find what he was looking for, he designed his own attachment and had it machined for this percussion device.

The Percussion Pro will arrive with the massage head separate and ready to attach (A). Push the attachment in firmly with the flat of your hand.

The Percussion Pro comes with a USB type-C charging cable that fits into any standard USB plug (B). The light on the Percussion Pro will turn red when it's charging, and blue when it's fully charged.

To turn on the Percussion Pro, hold the multi-function button for 2-3 seconds. Press again and the blue lights in the status indicator above the multi-function button will come on for each level as the oscillation increases. Press six times for the highest frequency. Using the device with 6 lights on produces an oscillation frequency of about 50 cps, which has been proven to reduce muscle tone and increase circulation. Press and hold the multi-function button for 2-3 seconds to turn off the Percussion Pro.



A



B

PLEASE NOTE: IF YOU HAVE A MEDICAL CONDITION OR HEALTH-RELATED PROBLEM, PLEASE CONSULT YOUR DOCTOR BEFORE USING THIS PRODUCT.

- DO NOT USE THIS PRODUCT ON OPEN WOUNDS, BLISTERS, OR BRUISES.
- DO NOT USE THIS PRODUCT IF YOU HAVE NEUROLOGIC CONDITIONS RESULTING IN LOSS OR ALTERED SENSATION
- DO NOT USE ON BONES, THIS PRODUCT IS FOR USE ON SOFT TISSUE ONLY
- DO NOT ALLOW USE OF PRODUCT BY MINORS: THIS IS NOT A TOY

Warning and Disclaimer

User assumes all risks associated with use of the PERCUSSION PRO. If you have any health problems listed above or you're unsure of whether to use this product, please consult your healthcare provider to determine whether use of this product is appropriate. Because injury may result from inappropriate use, such as the items in the warning list above, all of which is beyond control of the manufacturer and distributor, user assumes all risks. Under no circumstances shall the buyer be entitled to damages associated with the use of this product. Use of this product constitutes agreement to these terms.